



Measure Your Under-Bust:

- Measure around your ribcage, just under your bust, to get your band size. This measurement should be snug and level.
- Example: 95cm

Measure Bust Size:

- Measure from the centre of your chest (sternum) over the fullest part of your bust to your spine (non-operated breast), and then multiply this measurement by 2.
- Example: 56 inches × 2 = 112cm.
- Refer to the below chart to find out your Under bust and cup size

| EU underbust measurement in cm | 63-67 | 68-72 | 73-77 | 78-82 | 83-87 | 88-92 | 93-97 | 98-102 | 103-107 | 108-112 | 113-117 | 118-122 |
|-------------------------------------|-------|-------|-------|-------|-------|-------|-------|--------|---------|---------|---------|---------|
| L | + | + | + | + | + | + | + | + | + | + | + | ¥ |
| EU underbust circumference in cm | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 | 120 |

Sizes based on body measurements: 165–172 cm

2

Use the EU overbust measurement in cm to determine the corresponding cup size

| AA | 75-77 | 80-82 | 85-87 | 90-92 | 95-97 | 100-102 | 105-107 | 110-112 | 115-117 | 120-122 | 125-127 | 130-132 |
|----|-------|-------|--------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
| А | 77-79 | 82-84 | 87-89 | 92-94 | 97-99 | 102-104 | 107-109 | 112-114 | 117-119 | 122-124 | 127-129 | 132-134 |
| В | 79-81 | 84-86 | 89-91 | 94-96 | 99-101 | 104-106 | 109-111 | 114-116 | 119-121 | 124-126 | 129-131 | 134-136 |
| С | 81-83 | 86-88 | 91-93 | 96-98 | 101-103 | 106-108 | 111-113 | 116-118 | 121-123 | 126-128 | 131-133 | 136-138 |
| D | 83-85 | 88-90 | 93-95 | 98-100 | 103-105 | 108-110 | 113-115 | 118-120 | 123-125 | 128-130 | 133-135 | 138-140 |
| E | 85-87 | 90-92 | 95-97 | 100-102 | 105-107 | 110-112 | 115-117 | 120-122 | 125-127 | 130-132 | 135-137 | 140-142 |
| F | 87-89 | 92-94 | 97-99 | 102-104 | 107-109 | 112-114 | 117-119 | 122-124 | 127-129 | 132-134 | 137-139 | 142-144 |
| G | 89-91 | 94-96 | 99-101 | 104-106 | 109-111 | 114-116 | 119-121 | 124-126 | 129-131 | 134-136 | 139-141 | 144-146 |

• As per the above measurements your under bust size will be 95cm and cup size will be "C". As per the below chart the prosthesis size will be "9"

| Cup size | Underbust measurement | | | | | | | | | | | | | | | | | |
|----------|-----------------------|----|----|---------|----|----|----|--------|------------|------|-----|-----|-----|-------|-----|-----------------|-----|-----|
| AA | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 | 120 | | | | | | |
| А | 60 | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 | 120 | | | | | |
| В | | | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 | 120 | | | | |
| C | | | | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 | 120 | | | |
| D | | | | | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 | 120 | | |
| E | | | | | | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 | 120 | |
| F | | | | | | | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 | 120 |
| G | | | | | | | | | | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 |
| | | | | | | | | Overvi | ew of size | es 🖉 | | | | | | | | |
| | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | (9) | 10 | 11 | 12 | 13 | 14 | 15 ² | 16² | 172 |
| | | | | | | | | Prifo | rm sizes | | | | | | | | | |
| | | 1/ | 2 | 3/4 5/6 | | | | | 7/8 9/10 | | | 11/ | /12 | 13/14 | | | | |

HOW DO WE DETERMINE THE CORRECT BREAST FORM SIZE?



FX-A-2901, Marathon Futurex, Mafatlal Mills Compound, N. M. Joshi Marg, Lower Parel, Mumbai -400013 biosonixl@gmail.com | +91 8080901234